

# VEQ 360°



**Una risorsa indispensabile per confermare e migliorare la qualità di tutto il tuo laboratorio.**

Catalogo 2025



**IN STAND**

IN STAND International GMBH

**Programmi dedicati ad  
Allergologia, Autoimmunità e  
Proteine del complemento**



## Allergy diagnostics

## Allergy Diagnostic

Schellenberg / Zouboulis

720

Samples: 5 samples of 1.0 or 0.5 ml  
 Sample properties: Lyophilised serum  
 Shipping Information: Shipping at ambient temperature  
 Certificate Validity: 12 months

participation according to Rili-BÄK:  
 bi-annual

4 samples for determination specific IgE against inhaled allergens, food allergens and insect poison allergens

1 sample for determination of total IgE

	Jan. (1)	May (3)	July (4)	Oct. (6)
Registration deadline:	15.11.24	28.02.25	25.04.25	08.08.25
Sample Shipment:	14.01.25	29.04.25	24.06.25	07.10.25
Closing Date:	31.01.25	16.05.25	11.07.25	24.10.25

## Autoimmune diseases

## Autoimmune Diseases 01 - Connective Tissue Diseases (ANA)

Blüthner / Seelig

251

Samples: 2 samples of 0.3 ml  
 Sample properties: Plasma / Serum  
 Shipping Information: Shipping at ambient temperature  
 Certificate Validity: 12 months / 12 months

participation according to Rili-BÄK:  
 bi-annual

Antibodies against cell nuclei (antinuclear antibodies)

Autoantibodies against dsDNA  
 Collagenosis Screening

Autoantibodies against nuclei (ANA)  
 anti-DFS70 / LEDGF

	Mar. (2)	Aug. (5)
Registration deadline:	27.12.24	13.06.25
Sample Shipment:	25.02.25	12.08.25
Closing Date:	14.03.25	29.08.25

## Autoimmune Diseases 02 - Hepatic Syndromes

Blüthner / Seelig

253

Samples: 2 samples of 0.3 ml  
 Sample properties: Serum / Plasma  
 Shipping Information: Shipping at ambient temperature  
 Certificate Validity: 12 months

anti-Actin  
 anti-Mitochondria

anti-H+/K+-ATPase  
 anti-Parietal cells

anti-Microsomes (LKM 1)  
 anti-smooth muscle

	Mar. (2)	Aug. (5)
Registration deadline:	27.12.24	13.06.25
Sample Shipment:	25.02.25	12.08.25
Closing Date:	14.03.25	29.08.25

Samples: 2 samples of 0.3 ml Sample properties: Plasma / Serum Shipping Information: Shipping at ambient temperature Certificate Validity: 24 months		
Antibodies against ENA  anti-Jo 1 (Histidyl-tRNA-Synthetase) anti-SS-A / Ro anti-SS-A / Ro (60 kD - proteins) anti-Sm anti-Topoisomerase-1 (Scl 70)		
anti-PM / Scl anti-SS-A / Ro (52 kD - proteins) anti-SS-B / La anti-SmD anti-U1-snRNP		
		Aug. (5)
Registration deadline:		13.06.25
Sample Shipment:		12.08.25
Closing Date:		29.08.25

Samples: 2 samples of 0.3 ml Sample properties: Serum / Plasma Shipping Information: Shipping at ambient temperature Certificate Validity: 12 months / 12 months		participation according to Rili-BÄK: bi-annual
Antibodies against neutrophilic granulocytes (ANCA)  anti-MPO c-ANCA		
anti-PR3 p-ANCA	anti-glomerular basement membrane (GBM)	
	Mar. (2)	Aug. (5)
Registration deadline:	27.12.24	13.06.25
Sample Shipment:	25.02.25	12.08.25
Closing Date:	14.03.25	29.08.25

Samples: 2 samples of 0.3 ml Sample properties: Serum / Plasma Shipping Information: Shipping at ambient temperature Certificate Validity: 24 months			
anti- Desmosomes anti-Desmoglein 1 anti-Hemidesmosomes	anti-BPAg 180 anti-Desmoglein 3		anti-BPAg 230 anti-Desmoplakin I / II
Mar. (2)			
Registration deadline:	27.12.24		
Sample Shipment:	25.02.25		
Closing Date:	14.03.25		

Samples: 2 samples of 0.3 ml  
 Sample properties: Serum / Plasma  
 Shipping Information: Shipping at ambient temperature  
 Certificate Validity: 24 months

anti-Glutamic Acid Decarboxylase (GAD)	anti-Insulin (IAA)
anti-Islet cells (ICA)	anti-Islet cells - JDFU
anti-Tyrosine Phosphatase (IA2)	anti-Zinc Transporter 8 (ZnT8A)

**Aug. (5)**

Registration deadline:	13.06.25
Sample Shipment:	12.08.25
Closing Date:	29.08.25

Samples: 2 samples of 0.3 ml  
 Sample properties: Serum / Plasma  
 Shipping Information: Shipping at ambient temperature  
 Certificate Validity: 24 months

anti-P450 C17 - Hydroxylase	anti-P450 C21 - Hydroxylase
anti-P450 ssc - Hydroxylase	anti-adrenal cortex

**Aug. (5)**

Registration deadline:	13.06.25
Sample Shipment:	12.08.25
Closing Date:	29.08.25

Samples: 2 samples of 0.3 ml  
 Sample properties: Serum / Plasma  
 Shipping Information: Shipping at ambient temperature  
 Certificate Validity: 24 months

anti - Glia Cell Nuclear Antibody (AGNA)	anti-Amphiphysin 1
anti-HuD	anti-Neuronal nuclei
anti-Purkinje cells	anti-Ri / Nova 1
anti-SOX1	anti-Ta (Ma 2)
anti-Yo (CDR62)	

**Mar. (2)**

Registration deadline:	27.12.24
Sample Shipment:	25.02.25
Closing Date:	14.03.25

Samples: 2 samples of 0.3 ml Sample properties: Plasma / Serum Shipping Information: Shipping at ambient temperature Certificate Validity: 24 months		
anti-Asialoganglioside (AGA1) anti-Disialoganglioside (GD1b) anti-Monosialoganglioside (GM2) anti-Myelin-associated Glycoprotein (MAG) anti-Trisialoganglioside (GT1b)	anti-Disialoganglioside (GD1a) anti-Monosialoganglioside (GM1) anti-Monosialoganglioside (GM3) anti-Quadrosialoganglioside (GQ1b)	
Mar. (2)		
Registration deadline: Sample Shipment: Closing Date:	27.12.24 25.02.25 14.03.25	

## Autoimmune Diseases 10 - Myasthenia gravis

Samples: 2 samples of 0.3 ml Sample properties: Plasma / Serum Shipping Information: Shipping at ambient temperature Certificate Validity: 24 months		
anti-Acetylcholin receptors anti-Titin	anti-Musclespecific Receptor Tyrosinkinase (MuSK) anti-skeletal muscle	
Mar. (2)		
Registration deadline: Sample Shipment: Closing Date:	27.12.24 25.02.25 14.03.25	

## Autoimmune Diseases 11 - Gluten-Sensitive Enteropathies

Samples: 2 samples of 0.3 ml Sample properties: Plasma / Serum Shipping Information: Shipping at ambient temperature Certificate Validity: 12 months / 12 months			participation according to Rili-BÄK: bi-annual
anti-Endomysium IgA anti-Transglutaminase IgA	anti-Gliadin IgA anti-deamidated Gliadin IgA	anti-Gliadin IgG anti-deamidated Gliadin IgG	
<b>Mar. (2)</b>		<b>Aug. (5)</b>	
Registration deadline: Sample Shipment: Closing Date:	27.12.24 25.02.25 14.03.25	13.06.25 12.08.25 29.08.25	

Samples: 2 samples of 0.3 ml  
 Sample properties: Serum / Plasma  
 Shipping Information: Shipping at ambient temperature  
 Certificate Validity: 6 months / 12 months

participation according to Rili-BÄK:  
 quarterly

Rheumatoid Factor, Rheumatoid Factor (isotypes)  
 anti-MCV (mutated citrullinated Vimentin)

anti-CCP (cyclic citrulline peptide)

Additional terms see EQA No. 323!

	Mar. (2)	Aug. (5)
Registration deadline:	27.12.24	13.06.25
Sample Shipment:	25.02.25	12.08.25
Closing Date:	14.03.25	29.08.25

## Autoimmune Diseases 13 - Antiphospholipid Syndrome

Samples: 2 samples of 0.3 ml  
 Sample properties: Plasma / Serum  
 Shipping Information: Shipping at ambient temperature  
 Certificate Validity: 24 months

anti-Cardiolipin (IgA)      anti-Cardiolipin (IgG)      anti-Cardiolipin (IgM)  
 anti-β2-Glycoprotein (IgA)      anti-β2-Glycoprotein (IgG)      anti-β2-Glycoprotein (IgM)

	Aug. (5)
Registration deadline:	13.06.25
Sample Shipment:	12.08.25
Closing Date:	29.08.25

## Basophil Activation Test

Samples: 2 samples of 2.6 ml  
 Sample properties: EDTA whole blood, not stabilized, overnight shipping  
 Shipping Information: Cooled shipping  
 Certificate Validity: 12 months

Basophil activation by allergens bee and wasp. The stimulants are not included!

Due to high demand and limited sample material, we reserve the right to close orders early.

	May (3)	Oct. (6)
Registration deadline:	07.03.25	15.08.25
Sample Shipment:	06.05.25	14.10.25
Closing Date:	16.05.25	24.10.25

## Complement Analysis 01 – Auto-anti-C1q

Prohaszka / Kirschfink

245

Samples: 2 samples of 300 µl Sample properties: Serum Shipping Information: Shipping at ambient temperature Certificate Validity: 12 months		additional sample set not available	
Auto-anti C1q (IgG)	Auto-anti C1q (IgG) (qualitative)		
	<b>Mar. (2)</b>		<b>Nov. (6)</b>
Registration deadline:	27.12.24		08.08.25
Sample Shipment:	25.02.25	07.10.25	
Closing Date:	28.03.25	07.11.25	

## Complement Analysis 02 – Functions and proteins

Prohaszka / Kirschfink

246

Samples: 2 samples of 1 ml Sample properties: Lyophilised serum Shipping Information: Shipping at ambient temperature Certificate Validity: 12 months			Maximum order: 3 sets of samples
Alternative pathway C1q	C1 inhibitor: activity C3	C1 inhibitor: protein C4	
Classical pathway Lectin pathway	Factor H	Factor I	
Analysis of the total function of the 3 complement activation pathways (classical, alternative and lectin pathways) either by functional ELISA (WIELISA), by haemolytic (CH50, AH50) or by liposome-based assay (classical pathway). Participation with one parameter (e.g. CH50) is possible.			
	Mar. (2)	Nov. (6)	
Registration deadline:	27.12.24	08.08.25	
Sample Shipment:	25.02.25	07.10.25	
Closing Date:	28.03.25	07.11.25	

## Complement Analysis 03 – Complement activation products

Prohaszka / Kirschfink

247

Samples: 2 samples of 300 µl Sample properties: Lyophilised plasma Shipping Information: Shipping at ambient temperature Certificate Validity: 24 months	Maximum order: 3 sets of samples
C3a, C3d, Bb, sC5b-9	
Participation with one parameter (e.g. sC5b-9) possible	
	<b>Nov. (6)</b>
Registration deadline:	08.08.25
Sample Shipment:	07.10.25
Closing Date:	07.11.25

**Complement Analysis 04 – C3 – Nephritic Factor**

Prohaszka / Kirschfink

**248**

Samples: 2 samples of 300 µl  
Sample properties: Lyophilized serum  
Shipping Information: Shipping at ambient temperature  
Certificate Validity: 24 months

additional sample set not available

C3-nephritic factor (qualitative)

**Nov. (6)**

Registration deadline:	08.08.25
Sample Shipment:	07.10.25
Closing Date:	07.11.25

**Complement Analysis 05 – Auto-anti FH**

Prohaszka / Kirschfink

**249**

Samples: 2 samples of 300 µl  
Sample properties: Serum  
Shipping Information: Shipping at ambient temperature  
Certificate Validity: 24 months

additional sample set not available

Auto-anti FH-(qualitative)      Auto-anti-FH (quantitative)

**Nov. (6)**

Registration deadline:	08.08.25
Sample Shipment:	07.10.25
Closing Date:	07.11.25

**Complement Analysis 06 – Anti – C1 Inhibitor**

Prohaszka / Kirschfink

**250**

Samples: 2 samples of 300 µl  
Sample properties: Serum  
Shipping Information: Shipping at ambient temperature  
Certificate Validity: 24 months

additional sample set not available

Auto-anti-C1 inhibitor IgA (qualitative)      Auto-anti-C1 inhibitor IgG (qualitative)  
Auto-anti-C1 inhibitor IgM (qualitative)**Nov. (6)**

Registration deadline:	08.08.25
Sample Shipment:	07.10.25
Closing Date:	07.11.25