

VEQ 360°



Una risorsa indispensabile per confermare e migliorare la qualità di tutto il tuo laboratorio.

Catalogo 2025

LABQUALITY

Immunologia



Immunology

This program includes schemes for immunodiagnostic tests such as those for coeliac disease, rheumatoid factor and thyroid gland autoantibodies. For allergy diagnostics, review the allergology program in the clinical chemistry portfolio.

	1	2	3	4	5	6	7	8	9	10	11	12	
5935 ANCA and GbmAb Specimens: 2 liquid human serum or plasma samples, 0.5 mL. Examinations: Anti-neutrophil cytoplasmic Ab, Myeloperoxidase Ab, Proteinase-3 Ab and Glomerular basement membrane Ab.		●						●					EQA ³
3 Pre- and/or post-analytical cases in part of the rounds. Notes: Quantitative results are also processed (Pr3Ab, MPOAb, GbmAb).													
5900 Antinuclear antibodies Specimens: 3 liquid human serum or plasma samples, 0.6 mL. Examinations: ANA, ENAAb, RNPAAb, SmAb (SmDAB and/or SmBAB), SSAAb, SSBAb, Scl70Ab, CENP-B, CENP-A, Jo1Ab, dsDNA, HistAb, RibP Ab, RNApol				●						●			EQA ³
3 III Ab. Pre- and/or post-analytical cases in part of the rounds. Notes: Extractable antinuclear antigens and double-stranded deoxyribonucleic acid are included.													
5938 Autoimmune diagnostics, IFA interpretation Specimens: 3–5 cases (digital images).					●								
Examinations: Interpretation (ANA, ANCA and EMA images).													
5930 Autoimmune liver disease and gastric parietal cell antibodies Specimens: 2 liquid human serum or plasma samples, 0.4 mL.					●							●	
Examinations: Smooth muscle antibodies, Mitochondrial antibodies, Gastric parietal cell antibodies.													
5940 Coeliac disease, antibodies Specimens: 2 liquid human serum or plasma samples, 0.7 mL. Examinations: Endomysium antibodies, tissue transglutaminase antibodies, deamidated gliadin peptide antibodies, interpretation of the Total IgA concentration of the sample.		●				●					●		EQA ³
3 Pre- and/or post-analytical cases in part of the rounds. Notes: Quantitative results are also processed (tTGAbA, tTGAbG, DGPAbA, DGPAbG). Scheme is not suitable for POCT.													
5250 Interferon Gamma Release Assay (IGRA) for <i>Mycobacterium tuberculosis</i> Specimens: One sample set (contains 3 lyophilised samples, 1 liquid blank/ NIL sample and water to dissolve the samples) and a preanalytical case description including questions.		●			●			●				●	EQA ³
Examinations: Quantitative result and qualitative interpretation of TblNFg. The scheme is not suitable for the TB T-Spot test.													
5251 Interferon Gamma Release Assay (IGRA) for <i>Mycobacterium tuberculosis</i> - whole blood sample Specimens: 2 Li-heparin whole blood samples. Examinations: Quantitative result and qualitative interpretation of TblNFg.					●								NEW
3 Notes: EQA samples must be aliquoted into to test specific tubes (not provided) within 48h of phlebotomy. The samples are sent directly to the participating laboratories by the sample provider. By placing the order for this round, the participants give their consent to Labquality to transfer the participant contact information to the sample provider for shipping purposes.													
5937 Phospholipid antibodies Specimens: 2 liquid human serum or plasma samples, 0.5 mL. Examinations: Phospholipid antibodies, Cardiolipin antibodies (IgG and IgM), beta-2-glycoprotein antibodies (IgG and IgM).					●								
3 Notes: Quantitative results are also processed.													

1 2 3 4 5 6 7 8 9 10 11 12

5820 Rheumatoid factor and cyclic citrullinated peptide antibodies **3**

1	2	3	4	5	6	7	8	9	10	11	12
•			•			•			•		

Specimens: 2 liquid human-derived samples, 0.7 mL.

Examinations: Qualitative and quantitative RF, CCPAb.

1 2 3 4 5 6 7 8 9 10 11 12

5920 Thyroid gland antibodies **3**

1	2	3	4	5	6	7	8	9	10	11	12
		•			•				•		

Specimens: 2 liquid human serum or plasma samples, 0.4 mL.

Notes: Quantitative results are also processed.

Examinations: Thyroglobulin antibodies and thyroid peroxidase antibodies.
Pre- and/or post-analytical cases in part of the rounds.

1 2 3 4 5 6 7 8 9 10 11 12

5913 TSH receptor antibodies **3**

1	2	3	4	5	6	7	8	9	10	11	12
		•							•		

Specimens: 2 liquid human serum or plasma samples, 0.4 mL .

Notes: Quantitative results are also processed.

Examinations: Thyroid stimulating hormone receptor antibodies.

EGA³

LABQUALITY

**We are
Quality Makers**