

VEQ 360°



Una risorsa indispensabile per confermare e migliorare la qualità di tutto il tuo laboratorio.

Catalogo 2027

LABQUALITY

Anatomia patologica e Altre specialità



Clinical pathology

Seven high quality schemes are available for pathology laboratories. With changing topics in the rounds, both routine and more advanced needs are covered. The challenges are realistic and include also less commonly encountered clinically relevant cases. In the cytology and histopathology schemes, virtual microscopy is used. With this technology, viewing of several fields of vision and levels of focus are enabled on a computer screen simulating analysis with an optical microscope. Please see the new scheme for HPV-related head and neck cancer under Others >> Oncology.

Pathology, diagnostics

Clinical cytology - Gynaecological cytology (liquid based), virtual microscopy

Target: Laboratories

6701 Gynaecological cytology (liquid based), virtual microscopy												
Samples	Virtual images of at least 5 Papanicolaou stained slides of liquid based pap samples (ThinPrep). Diagnostics of cellular atypias in samples taken from gynaecological loci is assessed. Brief case histories and instructions are provided.											
Examinations	Observations and diagnoses.											
Rounds	JAN	FEB	MAR	APR	MAY	JUN	JUL	AUG	SEP	OCT	NOV	DEC
Accreditation	SFS-EN ISO/IEC 17043:2023, FINAS, PT02				Results/set		5					

Clinical cytology - Gynaecological cytology (smear), virtual microscopy

Target: Laboratories

6700 Gynaecological cytology (smear), virtual microscopy												
Samples	Virtual images of at least 5 Papanicolaou stained slides of conventional pap smear samples. The samples are selected from routine cytological material. Diagnostics of cellular atypias in samples taken from gynaecological loci is assessed. Brief case histories and instructions are provided.											
Examinations	Observations and diagnoses.											
Rounds	JAN	FEB	MAR	APR	MAY	JUN	JUL	AUG	SEP	OCT	NOV	DEC
Accreditation	SFS-EN ISO/IEC 17043:2023, FINAS, PT02				Results/set		5					

Clinical cytology - Non-gynaecological cytology, virtual microscopy

Target: Laboratories

6702 Non-gynaecological cytology, virtual microscopy												
Samples	Virtual images of Papanicolaou stained slides of nongynaecological cytosentrifuge or smear preparations or May-Grünwald- Giemsa stained smears. Images of at least 5 cases from representative loci. Brief case histories and instructions are provided.											
Examinations	Observations and diagnoses.											
Rounds	JAN	FEB	MAR	APR	MAY	JUN	JUL	AUG	SEP	OCT	NOV	DEC
Accreditation	SFS-EN ISO/IEC 17043:2023, FINAS, PT02					Results/set		5				

Histopathology, virtual microscopy

Target: Laboratories

6542 Histopathology, virtual microscopy												
Samples	Virtual images of at least 5 slides of miscellaneous tissue. Brief case histories and instructions are provided.											
Examinations	Topics: 1.2027 General Pathology 2.2027 Soft Tissue.											
Add. info	Topics may vary annually.											
Rounds	JAN	FEB	MAR	APR	MAY	JUN	JUL	AUG	SEP	OCT	NOV	DEC
Accreditation	SFS-EN ISO/IEC 17043:2023, FINAS, PT02					Results/set		5				

Pathology, technology

Histological staining techniques

Target: Laboratories

6543 Histological staining techniques												
Samples	Unstained paraffin embedded tissue from different tissue blocks or from one multiblock.											
Examinations	Topics 2027: 1/2027 Ziehl-Nielsen/AFB, PAS 2/2027 HE, Gram. Examinations: Staining of the slides. A set of stained slides is returned to Labquality for evaluation by an expert board.											
Add. info	Stains vary annually.											
Rounds	JAN	FEB	MAR	APR	MAY	JUN	JUL	AUG	SEP	OCT	NOV	DEC
Accreditation	SFS-EN ISO/IEC 17043:2023, FINAS, PT02					Results/set		1				

Immunohistochemical staining methods

Target: Laboratories

6600 Immunohistochemical staining methods												
Samples	Unstained paraffin embedded tissue from different tissue blocks or from one multiblock.											
Examinations	Topics 2027: 1/2027 PMS2, MLH1, MSH2, MSH6 2/2027 ER, PR, HER2 3/2027 BRAF V600E, Prame, Melan-A/MART-1. Examinations: Staining of the slides. A set of stained slides is returned to Labquality for evaluation by an expert board.											
Add. info	Stains vary annually.											
Rounds	JAN	FEB	MAR	APR	MAY	JUN	JUL	AUG	SEP	OCT	NOV	DEC
Accreditation	SFS-EN ISO/IEC 17043:2023, FINAS, PT02					Results/set		1				