

VEQ 360°



Una risorsa indispensabile per confermare e migliorare la qualità di tutto il tuo laboratorio.



Catalogo 2027

INSTAND

Complemento

EQAS/PT Properties Samples: 2 samples of 300 µl Sample Properties: Serum Max. number of sample sets: 1		EQAS/PT Experts: Prof. Dr. Zoltan Prohaszka Deputy EQAS/PT Experts: Prof. Dr. Michael Kirschfink			
Requirement for EQAS/PT according to Rili-BÄK: not applicable		Not subject to Rili-BÄK Validity Certificate: 24 months Parameter / Analytes: Auto-anti-C1 inhibitor IgA (qualitative) Auto-anti-C1 inhibitor IgG (qualitative) Auto-anti-C1 inhibitor IgM (qualitative)			
Shipping Information: Shipping not cooled		Safety Information: Process as patient sample material			
Dates	Q1		Q2	Q3	
					Nov.
Registration Deadline: Sample Shipment: Closing Date:					27-08-06 27-10-05 27-11-05



Complement Analysis - Anti-C1q Autoantibody

245

EQAS/PT Properties Samples: 2 samples of 300 µl Sample Properties: Serum Max. number of sample sets: 1			EQAS/PT Experts: Prof. Dr. Zoltan Prohaszka Deputy EQAS/PT Experts: Prof. Dr. Michael Kirschfink		
Requirement for EQAS/PT according to Rili-BÄK: not applicable			Not subject to Rili-BÄK Validity Certificate: 12 months Parameter / Analytes: Auto-anti C1q (IgG) Auto-anti C1q (IgG) (qualitative)		
Shipping Information: Shipping not cooled			Safety Information: Process as patient sample material		
Dates	Q1		Q2	Q3	
		Mar.			Nov.
Registration Deadline: Sample Shipment: Closing Date:		26-12-18 27-02-23 27-03-26			27-08-06 27-10-05 27-11-05

Complement Analysis - Anti-factor H (FH) Autoantibody

249

EQAS/PT Properties Samples: 2 samples of 300 µl Sample Properties: Serum Max. number of sample sets: 1			EQAS/PT Experts: Prof. Dr. Zoltan Prohaszka Deputy EQAS/PT Experts: Prof. Dr. Michael Kirschfink		
Requirement for EQAS/PT according to Rili-BÄK: not applicable			Not subject to Rili-BÄK Validity Certificate: 24 months Parameter / Analytes: Auto-anti FH-(qualitative) Auto-anti-FH (quantitative)		
Shipping Information: Shipping not cooled			Safety Information: Process as patient sample material		
Dates	Q1		Q2	Q3	
					Nov.
Registration Deadline: Sample Shipment: Closing Date:					27-08-06 27-10-05 27-11-05





Complement Analysis - C3 - Nephritic Factor

248

EQAS/PT Properties Samples: 2 samples of 300 µl Sample Properties: Lyophilised serum Max. number of sample sets: 1			EQAS/PT Experts: Prof. Dr. Zoltan Prohaszka Deputy EQAS/PT Experts: Prof. Dr. Michael Kirschfink		
Requirement for EQAS/PT according to Rili-BÄK: not applicable			Not subject to Rili-BÄK Validity Certificate: 24 months Parameter / Analytes: C3-nephritic factor (qualitative)		
Shipping Information: Shipping not cooled			Safety Information: Process as patient sample material		
Dates	Q1		Q2	Q3	
					Nov.
Registration Deadline: Sample Shipment: Closing Date:					27-08-06 27-10-05 27-11-05

Complement Analysis - Complement activation products

247

EQAS/PT Properties Samples: 2 samples of 300 µl Sample Properties: Lyophilised plasma Max. number of sample sets: 3 Additional Information: Participation with one parameter (e.g. sC5b-9) possible			EQAS/PT Experts: Prof. Dr. Zoltan Prohaszka Deputy EQAS/PT Experts: Prof. Dr. Michael Kirschfink		
Requirement for EQAS/PT according to Rili-BÄK: not applicable			Not subject to Rili-BÄK Validity Certificate: 24 months Parameter / Analytes: C3a, C3d, Bb, sC5b-9		
Shipping Information: Shipping not cooled			Safety Information: Process as patient sample material		
Dates	Q1		Q2	Q3	
					Nov.
Registration Deadline: Sample Shipment: Closing Date:					27-08-06 27-10-05 27-11-05



Complement Analysis - Functions (F) and Proteins (P)

246

EQAS/PT Properties Samples: 2 samples of 1 ml Sample Properties: Lyophilised serum Max. number of sample sets: 3 Additional Information: Analysis of the total function of the 3 complement activation pathways (classical, alternative and lectin pathways) either by functional ELISA (WIELISA), by haemolytic (CH50, AH50) or by liposome-based assay (classical pathway). Participation with one parameter (e.g. CH50) is possible.			EQAS/PT Experts: Prof. Dr. Zoltan Prohaszka Deputy EQAS/PT Experts: Prof. Dr. Michael Kirschfink		
Requirement for EQAS/PT according to Rili-BÄK: not applicable			Not subject to Rili-BÄK Validity Certificate: 12 months Parameter / Analytes: Classical pathway (F), Alternative pathway (F), Lectin pathway (FH), C1q, (P) C4 (P), C3 (P), C1 Inhibitor (F) (P), Factor H (P), Factor I (P)		
Shipping Information: Shipping not cooled			Safety Information: Process as patient sample material		
Dates	Q1		Q2	Q3	
		Mar.			Nov.
Registration Deadline:		26-12-18			27-08-06
Sample Shipment:		27-02-23			27-10-05
Closing Date:		27-03-26			27-11-05

