

VEQ 360°



Una risorsa indispensabile per confermare e migliorare la qualità di tutto il tuo laboratorio.



Catalogo 2027

INSTAND

Farmaci



Pharmaceuticals

Antiarrhythmics - Amiodaron

190

EQAS/PT Properties Samples: 2 samples of 3 ml Sample Properties: Lyophilised serum			EQAS/PT Experts: Prof. Dr. med. Werner Steimer Deputy EQAS/PT Experts: Dr. rer. nat. Heike Schneider		
Requirement for EQAS/PT according to Rili-BÄK: not applicable			Not subject to Rili-BÄK Validity Certificate: 12 months Parameter / Analytes: Amiodarone Desethylamiodarone		
Shipping Information: Shipping not cooled			Safety Information: Process as patient sample material		
Dates	Q1		Q2	Q3	
	Jan.		May	July	Oct.
Registration Deadline:	26-11-13		27-02-26	27-04-23	27-08-06
Sample Shipment:	27-01-12		27-04-27	27-06-22	27-10-05
Closing Date:	27-01-29		27-05-14	27-07-09	27-10-22

Antibacterials 01 - Teicoplanin

193

EQAS/PT Properties Samples: 2 samples of 2 ml Sample Properties: Lyophilised serum			EQAS/PT Experts: Prof. Dr. med. Werner Steimer Deputy EQAS/PT Experts: Dr. rer. nat. Heike Schneider		
Requirement for EQAS/PT according to Rili-BÄK: not applicable			Not subject to Rili-BÄK Validity Certificate: 12 months Parameter / Analytes: Teicoplanin		
Shipping Information: Shipping not cooled			Safety Information: Process as patient sample materia		
Dates	Q1		Q2	Q3	
	Jan.		May	July	Oct.
Registration Deadline:	26-11-13		27-02-26	27-04-23	27-08-06
Sample Shipment:	27-01-12		27-04-27	27-06-22	27-10-05
Closing Date:	27-01-29		27-05-14	27-07-09	27-10-22





Antibacterials 02

890

EQAS/PT Properties Samples: 2 samples of 1 ml Sample Properties: Lyophilised serum			EQAS/PT Experts: Dr. rer. nat. Heike Schneider Deputy EQAS/PT Experts: Prof. Dr. med. Werner Steimer		
Requirement for EQAS/PT according to Rili-BÄK: not applicable			Not subject to Rili-BÄK Validity Certificate: 12 months Parameter / Analytes: Ampicillin Meropenem Cefepime Piperacillin Ceftazidime Sulbactam Linezolid Tazobactam		
Shipping Information: Shipping not cooled			Safety Information: Process as patient sample material		
Dates	Q1		Q2	Q3	
	Jan.		May	Aug.	Oct.
Registration Deadline:	26-11-13		27-02-26	27-06-11	27-08-06
Sample Shipment:	27-01-12		27-04-27	27-08-10	27-10-05
Closing Date:	27-01-29		27-05-14	27-08-27	27-10-22

Anticonvulsants 01

191

EQAS/PT Properties Samples: 2 samples of 3 ml Sample Properties: Lyophilised serum			EQAS/PT Experts: Prof. Dr. med. Werner Steimer Deputy EQAS/PT Experts: Dr. rer. nat Heike Schneider		
Requirement for EQAS/PT according to Rili-BÄK: not applicable			Not subject to Rili-BÄK Validity Certificate: 12 months Parameter / Analytes: Lamotrigine Levetiracetam Sulthiame		
Shipping Information: Shipping not cooled			Safety Information: Process as patient sample material		
Dates	Q1		Q2	Q3	
	Jan.		May	Aug.	Oct.
Registration Deadline:	26-11-13		27-02-26	27-06-11	27-08-06
Sample Shipment:	27-01-12		27-04-27	27-08-10	27-10-05
Closing Date:	27-01-29		27-05-14	27-08-27	27-10-22



Anticonvulsants 02

192

EQAS/PT Properties Samples: 2 samples of 3 ml Sample Properties: Lyophilised serum			EQAS/PT Experts: Prof. Dr. med. Werner Steimer Deputy EQAS/PT Experts: Dr. rer. nat. Heike Schneider		
Requirement for EQAS/PT according to Rili-BÄK: not applicable			Not subject to Rili-BÄK Validity Certificate: 12 months Parameter / Analytes: 10-OH-Carbamazepine Carbamazepine-epoxide N-DesmethyImethsuximide		
Shipping Information: Shipping not cooled			Safety Information: Process as patient sample material		
Dates	Q1		Q2	Q3	
	Jan.		May		Aug.
Registration Deadline:	26-11-13		27-02-26		27-06-11
Sample Shipment:	27-01-12		27-04-27		27-08-10
Closing Date:	27-01-29		27-05-14		27-08-27
					27-10-05
					27-10-22

Anticonvulsants 03

871

EQAS/PT Properties Samples: 2 samples of 2.5 ml Sample Properties: Lyophilised serum			EQAS/PT Experts: Prof. Dr. med. Werner Steimer Deputy EQAS/PT Experts: Dr. rer. nat. Heike Schneider		
Requirement for EQAS/PT according to Rili-BÄK: not applicable			Not subject to Rili-BÄK Validity Certificate: 12 months Parameter / Analytes: Brivaracetam Eslicarbazepine Felbamate Lacosamide Perampanel Rufinamide		
Shipping Information: Shipping not cooled			Safety Information: Process as patient sample material		
Dates	Q1		Q2	Q3	
	Jan.		May		Aug.
Registration Deadline:	26-11-13		27-02-26		27-06-11
Sample Shipment:	27-01-12		27-04-27		27-08-10
Closing Date:	27-01-29		27-05-14		27-08-27
					27-10-05
					27-10-22





Anticonvulsants 04

872

EQAS/PT Properties Samples: 2 samples of 2.5 ml Sample Properties: Lyophilised serum			EQAS/PT Experts: Prof. Dr. med. Werner Steimer Deputy EQAS/PT Experts: Dr. rer. nat Heike Schneider		
Requirement for EQAS/PT according to Rili-BÄK: not applicable			Not subject to Rili-BÄK Validity Certificate: 12 months Parameter / Analytes: Gabapentin Pregabalin Tiagabine		
Shipping Information: Shipping not cooled			Safety Information: Process as patient sample material		
Dates	Q1		Q2	Q3	
	Jan.		May		Aug.
Registration Deadline:	26-11-13		27-02-26		27-06-11
Sample Shipment:	27-01-12		27-04-27		27-08-10
Closing Date:	27-01-29		27-05-14		27-08-27
					27-10-05
					27-10-22

Antifungals / Immunosuppressants

602

EQAS/PT Properties Samples: 2 samples of 1 mL Sample Properties: Lyophilised plasma			EQAS/PT Experts: Prof. Dr. med. Nicolas von Ahsen Deputy EQAS/PT Experts: Dr. Patricia Kaiser (PhD)		
Requirement for EQAS/PT according to Rili-BÄK: not applicable			Not subject to Rili-BÄK Validity Certificate: 12 months Parameter / Analytes: Hydroxy-Itraconazole Itraconazole Mycophenolic acid		
Shipping Information: Shipping not cooled			Safety Information: Process as patient sample material		
Dates	Q1		Q2	Q3	
		Mar.			Oct.
Registration Deadline:		26-12-18			27-08-06
Sample Shipment:		27-02-23			27-10-05
Closing Date:		27-03-12			27-10-22



Antifungals and Chloramphenicol

198

EQAS/PT Properties Samples: 2 samples of 5 ml (+ 2 samples of 3 ml antibiotics) Sample Properties: Lyophilised serum			EQAS/PT Experts: Prof. Dr. med. Werner Steimer Deputy EQAS/PT Experts: Dr. rer. nat Heike Schneider		
Requirement for EQAS/PT according to Rili-BÄK: not applicable			Not subject to Rili-BÄK Validity Certificate: 12 months Parameter / Analytes: Chloramphenicol Posaconazole Fluconazole Voriconazole Flucytosine		
Shipping Information: Shipping not cooled			Safety Information: Process as patient sample material		
Dates	Q1		Q2	Q3	
			May		Oct.
Registration Deadline: Sample Shipment: Closing Date:			27-02-26 27-04-27 27-05-14		27-08-06 27-10-05 27-10-22

Antipsychotics/Antidepressives 01

200

EQAS/PT Properties Samples: 2 samples of 2 ml Sample Properties: Lyophilised serum			EQAS/PT Experts: Prof. Dr. med. Werner Steimer Deputy EQAS/PT Experts: Dr. rer. nat Heike Schneider		
Requirement for EQAS/PT according to Rili-BÄK: not applicable			Not subject to Rili-BÄK Validity Certificate: 12 months Parameter / Analytes: Clozapine N-Desmethylozapine		
Shipping Information: Shipping not cooled			Safety Information: Process as patient sample material		
Dates	Q1		Q2	Q3	
	Jan.	Mar.	May	July	Aug.
Registration Deadline: Sample Shipment: Closing Date:	26-11-13 27-01-12 27-01-29	26-12-18 27-02-23 27-03-12	27-02-26 27-04-27 27-05-14	27-04-23 27-06-22 27-07-09	27-06-11 27-08-10 27-08-27





Antipsychotics/Antidepressives 02

873

EQAS/PT Properties Samples: 2 samples of 2.5 ml Sample Properties: Lyophilised serum				EQAS/PT Experts: Prof. Dr. med. Werner Steimer Deputy EQAS/PT Experts: Dr. rer. nat Heike Schneider		
Requirement for EQAS/PT according to Rili-BÄK: not applicable				Not subject to Rili-BÄK Validity Certificate: 12 months Parameter / Analytes: Amisulpride		

Antipsychotics/Antidepressives 03

874

EQAS/PT Properties Samples: 2 samples of 2.5 ml Sample Properties: Lyophilised serum				EQAS/PT Experts: Prof. Dr. med. Werner Steimer Deputy EQAS/PT Experts: Dr. rer. nat Heike Schneider		
Requirement for EQAS/PT according to Rili-BÄK: not applicable				Not subject to Rili-BÄK Validity Certificate: 12 months Parameter / Analytes: Desmethylfluoxetine Mianserine Fluoxetine Perazine		
Shipping Information: Shipping not cooled				Safety Information: Process as patient sample material		
Dates	Q1		Q2	Q3		Q4
	Jan.	Mar.	May	July	Aug.	Oct.
Registration Deadline:	26-11-13	26-12-18	27-02-26	27-04-23	27-06-11	27-08-06
Sample Shipment:	27-01-12	27-02-23	27-04-27	27-06-22	27-08-10	27-10-05
Closing Date:	27-01-29	27-03-12	27-05-14	27-07-09	27-08-27	27-10-22



Antipsychotics/Antidepressives 04

875

EQAS/PT Properties Samples: 2 samples of 2.5 ml Sample Properties: Lyophilised serum				EQAS/PT Experts: Prof. Dr. med. Werner Steimer Deputy EQAS/PT Experts: Dr. rer. nat Heike Schneider		
Requirement for EQAS/PT according to Rili-BÄK: not applicable				Not subject to Rili-BÄK Validity Certificate: 12 months Parameter / Analytes: Cariprazine		

Antivirals

850

EQAS/PT Properties Samples: 2 samples of 1 - 3 ml Sample Properties: Lyophilised serum Additional Information: Please note that some analytes may interfere with the HPLC methods with UV detection. For this reason, the schedule-based composition will be adjusted as shown in the following pdf-download-link: https://www.instand-ev.de/EQA850-Analytes				EQAS/PT Experts: Dr. rer. nat Heike Schneider Deputy EQAS/PT Experts: Prof. Dr. med. Werner Steimer																						
Requirement for EQAS/PT according to Rili-BÄK: not applicable			Not subject to Rili-BÄK Validity Certificate: 12 months Parameter / Analytes: <table><tr><td>Aciclovir</td><td>Efavirenz</td><td>Raltegravir</td></tr><tr><td>Atazanavir</td><td>Elvitegravir</td><td>Ribavirin</td></tr><tr><td>Bictegravir</td><td>Etravirine</td><td>Rilpivirin</td></tr><tr><td>Cabotegravir</td><td>Ganciclovir</td><td>Ritonavir T</td></tr><tr><td>Darunavir</td><td>Lopinavir</td><td>enofovir</td></tr><tr><td>Dolutegravir</td><td>Maraviroc</td><td></td></tr><tr><td>Doravirine</td><td>Nevirapine</td><td></td></tr></table>			Aciclovir	Efavirenz	Raltegravir	Atazanavir	Elvitegravir	Ribavirin	Bictegravir	Etravirine	Rilpivirin	Cabotegravir	Ganciclovir	Ritonavir T	Darunavir	Lopinavir	enofovir	Dolutegravir	Maraviroc		Doravirine	Nevirapine	
Aciclovir	Efavirenz	Raltegravir																								
Atazanavir	Elvitegravir	Ribavirin																								
Bictegravir	Etravirine	Rilpivirin																								
Cabotegravir	Ganciclovir	Ritonavir T																								
Darunavir	Lopinavir	enofovir																								
Dolutegravir	Maraviroc																									
Doravirine	Nevirapine																									
Shipping Information: Shipping not cooled			Safety Information: Process as patient sample material																							
Dates	Q1		Q2	Q3		Q4																				
		Mar.	May		Aug.	Oct.																				
Registration Deadline:		26-12-18	27-02-26		27-06-11	27-08-06																				
Sample Shipment:		27-02-23	27-04-27		27-08-10	27-10-05																				
Closing Date:		27-03-12	27-05-14		27-08-27	27-10-22																				



876

EQAS/PT Properties Samples: 2 samples of 2.5 ml Sample Properties: Lyophilised serum				EQAS/PT Experts: Prof. Dr. med. Werner Steimer Deputy EQAS/PT Experts: Dr. rer. nat Heike Schneider			
Requirement for EQAS/PT according to Rili-BÄK: not applicable				Not subject to Rili-BÄK Validity Certificate: 12 months Parameter / Analytes: Alprazolam Chlordiazepoxide Desalkylflurazepam Diazepam Nordiazepam Oxazepam Trazodone			
Shipping Information: Shipping not cooled				Safety Information: Process as patient sample material			
Dates	Q1		Q2	Q3		Q4	
	Jan.		May		Aug.	Oct.	
Registration Deadline:	26-11-13		27-02-26		27-06-11	27-08-06	
Sample Shipment:	27-01-12		27-04-27		27-08-10	27-10-05	
Closing Date:	27-01-29		27-05-14		27-08-27	27-10-22	

877

EQAS/PT Properties Samples: 2 samples of 2.5 ml Sample Properties: Lyophilised serum				EQAS/PT Experts: Prof. Dr. med. Werner Steimer Deputy EQAS/PT Experts: Dr. rer. nat Heike Schneider		
Requirement for EQAS/PT according to Rili-BÄK: not applicable				Not subject to Rili-BÄK Validity Certificate: 12 months Parameter / Analytes: Clonazepam Nitrazepam Flunitrazepam Temazepam Flurazepam		
Shipping Information: Shipping not cooled				Safety Information: Process as patient sample material		
Dates	Q1		Q2	Q3		Q4
	Jan.		May		Aug.	Oct.
Registration Deadline: Sample Shipment: Closing Date:	26-11-13 27-01-12 27-01-29		27-02-26 27-04-27 27-05-14		27-06-11 27-08-10 27-08-27	27-08-06 27-10-05 27-10-22



Benzodiazepines 03

878

EQAS/PT Properties Samples: 2 samples of 2.5 ml Sample Properties: Lyophilised serum			EQAS/PT Experts: Prof. Dr. med. Werner Steimer Deputy EQAS/PT Experts: Dr. rer. nat Heike Schneider			
Requirement for EQAS/PT according to Rili-BÄK: not applicable			Not subject to Rili-BÄK Validity Certificate: 12 months Parameter / Analytes: Bromazepam Clobazam (Clobazepam) Lorazepam			
Shipping Information: Shipping not cooled			Safety Information: Process as patient sample material			
Dates	Q1		Q2	Q3		Q4
	Jan.		May		Aug.	Oct.
Registration Deadline:	26-11-13		27-02-26		27-06-11	27-08-06
Sample Shipment:	27-01-12		27-04-27		27-08-10	27-10-05
Closing Date:	27-01-29		27-05-14		27-08-27	27-10-22

Immunosuppressants

601

EQAS/PT Properties Samples: 4 samples of 1 ml Sample Properties: Lyophilised whole blood			EQAS/PT Experts: Prof. Dr. med. Nicolas von Ahsen Deputy EQAS/PT Experts: Dr. Patricia Kaiser (PhD)		
Requirement for EQAS/PT according to Rili-BÄK: not applicable			Not subject to Rili-BÄK Validity Certificate: 12 months Parameter / Analytes: Cyclosporine A		





Paracetamol and NSAIDs

194

EQAS/PT Properties Samples: 4 samples of 1 ml Sample Properties: Lyophilised serum			EQAS/PT Experts: Prof. Dr. med. Werner Steimer Deputy EQAS/PT Experts: Dr. rer. nat Heike Schneider			
Requirement for EQAS/PT according to Rili-BÄK: not applicable			Not subject to Rili-BÄK Validity Certificate: 12 months Parameter / Analytes: Paracetamol Salicylic acid			
Shipping Information: Shipping not cooled			Safety Information: Process as patient sample material			
Dates	Q1		Q2	Q3		Q4
	Jan.		May	July		Oct.
Registration Deadline:	26-11-13		27-02-26	27-04-23		27-08-06
Sample Shipment:	27-01-12		27-04-27	27-06-22		27-10-05
Closing Date:	27-01-29		27-05-14	27-07-09		27-10-22

Psychostimulants - ADHS-Treatment

879

EQAS/PT Properties Samples: 2 samples of 2.5 ml Sample Properties: Lyophilised serum			EQAS/PT Experts: Prof. Dr. med. Werner Steimer Deputy EQAS/PT Experts: Dr. rer. nat Heike Schneider			
Requirement for EQAS/PT according to Rili-BÄK: not applicable			Not subject to Rili-BÄK Validity Certificate: 12 months Parameter / Analytes: Atomoxetine Methylphenidate Ritalinic acid			
Shipping Information: Shipping not cooled			Safety Information: Process as patient sample material			
Dates	Q1		Q2	Q3		Q4
	Jan.		May		Aug.	Oct.
Registration Deadline:	26-11-13		27-02-26		27-06-11	27-08-06
Sample Shipment:	27-01-12		27-04-27		27-08-10	27-10-05
Closing Date:	27-01-29		27-05-14		27-08-27	27-10-22



Psychotherapeutic Drugs 01

860

EQAS/PT Properties Samples: 2 samples of 5 ml Sample Properties: Lyophilised serum			EQAS/PT Experts: Prof. Dr. med. Dr. rer. nat. Dipl. Biochem. Ekkehard Haen Deputy EQAS/PT Experts: Florian Ridders			
Requirement for EQAS/PT according to Rili-BÄK: not applicable			Not subject to Rili-BÄK Validity Certificate: 12 months Parameter / Analytes: 10-OH-Carbamazepine O-Desmethylvenlafaxine Aripiprazole Perazine Citalopram Quetiapine Lamotrigine Venlafaxine			
Shipping Information: Shipping not cooled			Safety Information: Process as patient sample material			
Dates	Q1		Q2	Q3		Q4
	Jan.		May	July		Oct.
Registration Deadline:	26-11-13		27-02-26	27-04-23		27-08-06
Sample Shipment:	27-01-12		27-04-27	27-06-22		27-10-05
Closing Date:	27-01-29		27-05-14	27-07-09		27-10-22

Psychotherapeutic Drugs 02

861

EQAS/PT Properties Samples: 2 samples of 5 ml Sample Properties: Lyophilised serum			EQAS/PT Experts: Prof. Dr. med. Dr. rer. nat. Dipl. Biochem. Ekkehard Haen Deputy EQAS/PT Experts: Florian Ridders			
Requirement for EQAS/PT according to Rili-BÄK: not applicable			Not subject to Rili-BÄK Validity Certificate: 12 months Parameter / Analytes: Amisulpride Escitalopram Mirtazapine O-Desmethylvenlafaxine Quetiapine Venlafaxine			
Shipping Information: Shipping not cooled			Safety Information: Process as patient sample material			
Dates	Q1		Q2	Q3		Q4
	Jan.		May	July		Oct.
Registration Deadline:	26-11-13		27-02-26	27-04-23		27-08-06
Sample Shipment:	27-01-12		27-04-27	27-06-22		27-10-05
Closing Date:	27-01-29		27-05-14	27-07-09		27-10-22





Psychotherapeutic Drugs 03

862

EQAS/PT Properties Samples: 2 samples of 5 ml Sample Properties: Lyophilised serum			EQAS/PT Experts: Prof. Dr. med. Dr. rer. nat. Dipl. Biochem. Ekkehard Haen Deputy EQAS/PT Experts: Florian Ridders		
Requirement for EQAS/PT according to Rili-BÄK: not applicable			Not subject to Rili-BÄK Validity Certificate: 12 months Parameter / Analytes: 9-OH-Risperidone / Paliperidone Olanzapine Clozapine Quetiapine Duloxetine Risperidone Haloperidol Ziprasidone N-Desmethylozapine		
Shipping Information: Shipping not cooled			Safety Information: Process as patient sample material		
Dates	Q1		Q2	Q3	
	Jan.		May	July	Oct.
Registration Deadline:	26-11-13		27-02-26	27-04-23	27-08-06
Sample Shipment:	27-01-12		27-04-27	27-06-22	27-10-05
Closing Date:	27-01-29		27-05-14	27-07-09	27-10-22

Psychotherapeutic Drugs 04

863

EQAS/PT Properties Samples: 2 samples of 5 ml Sample Properties: Lyophilised serum			EQAS/PT Experts: Prof. Dr. med. Dr. rer. nat. Dipl. Biochem. Ekkehard Haen Deputy EQAS/PT Experts: Florian Ridders		
Requirement for EQAS/PT according to Rili-BÄK: not applicable			Not subject to Rili-BÄK Validity Certificate: 12 months Parameter / Analytes: Amitriptyline Imipramine Clomipramine Maprotiline Clozapine N-Desmethylozapine Desipramine Norclomipramine Doxepin Nordoxepin		
Shipping Information: Shipping not cooled			Safety Information: Process as patient sample material		
Dates	Q1		Q2	Q3	
	Jan.		May	July	Oct.
Registration Deadline:	26-11-13		27-02-26	27-04-23	27-08-06
Sample Shipment:	27-01-12		27-04-27	27-06-22	27-10-05
Closing Date:	27-01-29		27-05-14	27-07-09	27-10-22



Respiratory Stimulants - Caffeine

870

EQAS/PT Properties Samples: 2 samples of 2.5 ml Sample Properties: Lyophilised serum			EQAS/PT Experts: Prof. Dr. med. Werner Steimer Deputy EQAS/PT Experts: Dr. rer. nat Heike Schneider		
Requirement for EQAS/PT according to Rili-BÄK: not applicable			Not subject to Rili-BÄK Validity Certificate: 12 months Parameter / Analytes: Caffeine		
Shipping Information: Shipping not cooled			Safety Information: Process as patient sample material		
Dates	Q1		Q2	Q3	
			May		Oct.
Registration Deadline: Sample Shipment: Closing Date:			27-02-26 27-04-27 27-05-14		27-08-06 27-10-05 27-10-22

Therapeutic drugs - Antibiotics, Digitoxin, MTX

197

EQAS/PT Properties Samples: 2 samples of 5 ml Sample Properties: Lyophilised serum			EQAS/PT Experts: Prof. Dr. med. Werner Steimer Deputy EQAS/PT Experts: Dr. rer. nat Heike Schneider Dr. Patricia Kaiser (PhD)		
Requirement for EQAS/PT according to Rili-BÄK: quarterly Validity Certificate: 6 months Parameter / Analytes: Digitoxin Vancomycin			Not subject to Rili-BÄK Validity Certificate: 12 months Parameter / Analytes: Amicacin Methotrexate Gentamicin Tobramycin		
Shipping Information: Shipping not cooled			Safety Information: Process as patient sample material		
Dates	Q1		Q2	Q3	
	Jan.	Mar.	May	July	Aug.
Registration Deadline: Sample Shipment: Closing Date:	26-11-13 27-01-12 27-01-29	26-12-18 27-02-23 27-03-12	27-02-26 27-04-27 27-05-14	27-04-23 27-06-22 27-07-09	27-06-11 27-08-10 27-08-27





Therapeutic drugs - Antiepileptics, Digoxin, Lithium and others

195

EQAS/PT Properties Samples: 2 samples of 5 ml Sample Properties: Lyophilised serum				EQAS/PT Experts: Prof. Dr. med. Werner Steimer Deputy EQAS/PT Experts: Dr. rer. nat Heike Schneider Dr. Patricia Kaiser (PhD)		
Requirement for EQAS/PT according to Rili-BÄK: quarterly Validity Certificate: 6 months Parameter / Analytes: Carbamacepine Phenytoin Lithium Theophylline Phenobarbital Valproic Acid				Not subject to Rili-BÄK Validity Certificate: 12 months Parameter / Analytes: Digoxin Ethosuximid Primidon		
Shipping Information: Shipping not cooled				Safety Information: Process as patient sample material		
Dates	Q1		Q2	Q3		Q4
	Jan.	Mar.	May	July	Aug.	Oct.
Registration Deadline:	26-11-13	26-12-18	27-02-26	27-04-23	27-06-11	27-08-06
Sample Shipment:	27-01-12	27-02-23	27-04-27	27-06-22	27-08-10	27-10-05
Closing Date:	27-01-29	27-03-12	27-05-14	27-07-09	27-08-27	27-10-22

Tri-/Tetracyclic Antidepressants 03 - Quantitative in Serum

202

EQAS/PT Properties Samples: 2 samples of 2.5 ml Sample Properties: Lyophilised serum			EQAS/PT Experts: Prof. Dr. med. Werner Steimer Deputy EQAS/PT Experts: Dr. rer. nat Heike Schneider		
Requirement for EQAS/PT according to Rili-BÄK: not applicable			Not subject to Rili-BÄK Validity Certificate: 12 months Parameter / Analytes: Desipramine		



Tricyclic Antidepressants 01 - Quantitative in Serum

199

EQAS/PT Properties Samples: 2 samples of 2.5 ml Sample Properties: Lyophilised serum			EQAS/PT Experts: Prof. Dr. med. Werner Steimer Deputy EQAS/PT Experts: Dr. rer. nat Heike Schneider		
Requirement for EQAS/PT according to Rili-BÄK: not applicable			Not subject to Rili-BÄK Validity Certificate: 12 months Parameter / Analytes: Amitriptyline Norclomipramine Clomipramine Nortriptyline Tricyclic Antidepressants - Screen		
Shipping Information: Shipping not cooled			Safety Information: Process as patient sample material		
Dates	Q1		Q2	Q3	
	Jan.		May	July	Oct.
Registration Deadline:	26-11-13		27-02-26	27-04-23	27-08-06
Sample Shipment:	27-01-12		27-04-27	27-06-22	27-10-05
Closing Date:	27-01-29		27-05-14	27-07-09	27-10-22

Tricyclic Antidepressants 02 - Quantitative in Serum

201

EQAS/PT Properties Samples: 2 samples of 2.5 ml Sample Properties: Lyophilised serum			EQAS/PT Experts: Prof. Dr. med. Werner Steimer Deputy EQAS/PT Experts: Dr. rer. nat Heike Schneider		
Requirement for EQAS/PT according to Rili-BÄK: not applicable			Not subject to Rili-BÄK Validity Certificate: 12 months Parameter / Analytes: Doxepin Nordoxepin Nortrimipramine Tricyclic Antidepressants - Screen Trimipramine		
Shipping Information: Shipping not cooled			Safety Information: Process as patient sample material		
Dates	Q1		Q2	Q3	
	Jan.		May	July	Oct.
Registration Deadline:	26-11-13		27-02-26	27-04-23	27-08-06
Sample Shipment:	27-01-12		27-04-27	27-06-22	27-10-05
Closing Date:	27-01-29		27-05-14	27-07-09	27-10-22





Tricyclic Antidepressants 04 - Qualitative in Serum

203

EQAS/PT Properties**Samples:** 2 samples of 2.5 ml**Sample Properties:** Lyophilised serum**Additional Information:**

- May contain the following substances: Amitriptyline, Clomipramine, Desipramine, Doxepin, Imipramine, Norclomipramine, Nordoxepin, Nortrimipramine, Nortriptyline, Trimipramine
- For POCT systems: Please ensure that the sample matrix used in the EQA scheme is suitable for your test system. In case of uncertainty, we recommend contacting the manufacturer or supplier of your device.

**EQAS/PT Experts:**

Prof. Dr. med. Werner Steimer

Deputy EQAS/PT Experts:

Dr. rer. nat Heike Schneider

Requirement for EQAS/PT according to Rili-BÄK: quarterly**Validity Certificate:** 6 months**Parameter / Analytes:**

Tricyclic Antidepressants (qualitative)- Screening

Not subject to Rili-BÄK:**Validity Certificate:** 12 months**Parameter / Analytes:**

Amitriptyline	Norclomipramine
Clomipramine	Nordoxepin
Desipramine	Nortrimipramine
Doxepin	Nortriptyline
Imipramine	Trimipramine

Shipping Information:

Shipping not cooled

Safety Information:

Process as patient sample material

Dates	Q1		Q2	Q3		Q4
	Jan.	Mar.	May	July	Aug.	Oct.
Registration Deadline:	26-11-13	26-12-18	27-02-26	27-04-23	27-06-11	27-08-06
Sample Shipment:	27-01-12	27-02-23	27-04-27	27-06-22	27-08-10	27-10-05
Closing Date:	27-01-29	27-03-12	27-05-14	27-07-09	27-08-27	27-10-22