

Lenti RapidTiter™ Sticks


Designed for semi-quantitative detection of the lentiviral P24 protein, Lenti RapidTiter™ Sticks revolutionize the way researchers measure lentiviral titer. Accurate titer determination allows researchers to fine-tune transduction conditions, ensuring optimal viral concentrations for efficient gene delivery.

Key Feature and Benefits

- **Fast Results:** Load sample, wait, and compare the band intensity, in just 10 minutes!
- **Broad Detection Range:** Detect titers between 2.5×10^5 - 2.5×10^7 TU/mL (*Note: higher titers can be assayed with dilution*)
- **User-Friendly Design:** Simple design and easy protocol allowing for seamless testing for researchers at various skill levels.

Name	Size	Catalog #
Lenti RapidTiter™ Sticks	5 tests	EB200001
	20 tests	EB200002
	50 tests	EB200003

Estimate your Titer in 4 Easy Steps



Browse related
lentiviral products



CONTACT US:

P 301.340.3188

F 301.340.9254

E sales@origene.com