



JusChek⁺

PRODUCT CATALOGUE



ACRO BIOTECH INC.
www.acrobiotech.com

ABOUT US & CERTIFICATES

Acro Biotech is a California-based manufacturer of rapid tests offering the products from state-of-the-art ISO13485:2016 accredited CMDCAS-compliant manufacturing facilities. Our core team comprises of seasoned professionals in the field of rapid tests with more than two decades of experience in reputed rapid test manufacturing companies, working in senior positions. Our product range includes USFDA (510)K approved DOA tests as well as an expansive range comprising of women's health related rapid tests, cardiac markers, tumor markers, infectious disease and veterinary tests for international market.


Our tests work with a wide range of specimens, including urine, feces, oral fluid, whole blood, serum, plasma, hair and swabs.

At the cusp of innovation and commercialization, we are working hard to introduce new products in the market that would set the future growth path, including the combination of digitization and ease of tests with lateral flow readers and fluorescence immunoassay analyzers.

EU ISO 13485




US FDA



EU CE1434 OTC



Health Canada



UK MHRA



Saudi Arabia MDMA



Manufacturing License (GMP)



CONTENTS

Rapid Test

Women's Health	-----	1-2
Infectious Disease	-----	3-10
Drugs of Abuse	-----	11-16
Tumor Markers	-----	17
Cardiac Markers	-----	17
Others	-----	18



Rapid Tests for Women's Health

Pregnancy

Pregnancy (hCG)

Specimen: Urine, WB/S/P

Format: Dipstick, Cassette, Midstream, Digital Midstream

Cut-Off: 10 mIU/mL, 20 mIU/mL, 25 mIU/mL*, 40 mIU/mL, 200 mIU/mL, 500 mIU/mL, 2000 mIU/mL

Pack: 1 T, 25 T, 50 T

CE Status: CE, US FDA



Ovulation

AMH (Anti-Mullerian Hormone)

Specimen: WB/S/P

Format: Cassette

Cut-Off: 1-7 ng/mL

Pack: 10 T

CE Status: CE



Ovulation (LH)

Specimen: Urine

Format: Dipstick, Cassette, Midstream

Cut-Off: 30 mIU/mL, 35 mIU/mL, 40 mIU/mL*

Pack: 5 T, 25 T, 50 T

CE Status: CE, US FDA

Menopause

FSH (Follicle Stimulating Hormone)

Specimen: Urine

Format: Dipstick, Cassette, Midstream

Cut-Off: 25 mIU/mL

Pack: 2 T, 25 T, 50 T

CE Status: CE

* Standard Cut-off

Preterm Delivery (PROM)

Fetal Fibronectin (fFN)

Specimen: Vaginal Secretion

Format: Cassette

Pack: 20 T

CE Status: CE

Insulin-like Growth Factor-binding Protein 1 (IGFBP-1)

Specimen: Vaginal Secretion

Format: Cassette

Pack: 20 T

CE Status: CE



ToRCH Range

Toxo IgG/IgM

Specimen: WB/S/P

Format: Cassette

Pack: 25 T

CE Status: Non-CE

Rubella IgG/IgM

Specimen: WB/S/P

Format: Cassette

Pack: 25 T

CE Status: Non-CE

CMV IgG/IgM

Specimen: WB/S/P

Format: Cassette

Pack: 25 T

CE Status: Non-CE

HSV 1/2 IgG/IgM

Specimen: WB/S/P

Format: Cassette

Pack: 25 T

CE Status: CE



ToRCH IgG/IgM

Specimen: WB/S/P

Format: Cassette


Pack: 25 T

CE Status: Non-CE

Rapid Tests for Infectious Diseases


Gastrointestinal Diseases

Cat No.	Product	Specimen	Format	Pack	CE Status
IAD-602	Adenovirus	Feces	Cassette	25 T	CE
IRO-602	Rotavirus	Feces	Cassette	25 T	CE
INO-602	Norovirus	Feces	Cassette	10 T	CE
IRAC-625	Rotavirus and Adenovirus Combo	Feces	Cassette	25 T	CE
INRA-635	Norovirus, Rotavirus and Adenovirus Combo	Feces	Cassette	10 T	CE
IAS-602	Astrovirus	Feces	Cassette	10 T	CE
IMVD-645	Norovirus, Rotavirus, Adenovirus and Astrovirus Combo	Feces	Cassette	10 T	CE
IEH-602	<i>Entamoeba histolytica</i>	Feces	Cassette	10 T	CE
IGL-602	<i>Giardia Lamblia</i>	Feces	Cassette	10 T	CE
ICR-602	<i>Cryptosporidium</i>	Feces	Cassette	10 T	CE
ICGC-625	<i>Cryptosporidium</i> and <i>Giardia Lamblia</i> Combo	Feces	Cassette	10 T	CE
IEGC-602	<i>Entamoeba/Giardia/Crypto</i> (1 Window)	Feces	Cassette	10 T	CE
IEGC-635	<i>Entamoeba/Giardia/Crypto</i> Combo (3 Windows)	Feces	Cassette	10 T	CE
ICAM-602	<i>Campylobacter</i>	Feces	Cassette	10 T	CE
ICDG-602	<i>Clostridium difficile</i> GDH	Feces	Cassette	10 T	CE
ICDT-625B	<i>Clostridium difficile</i> Toxin A +Toxin B Combo	Feces	Cassette	10 T	CE
ICD-635B	<i>Clostridium difficile</i> GDH+ Toxin A +Toxin B Combo	Feces	Cassette	25 T	CE
IHP-302	<i>H. pylori</i> Antibody	S/P	Cassette	25 T	CE
IHP-401	<i>H. pylori</i> Antibody	WB/S/P	Dipstick	50 T	CE
IHP-402	<i>H. pylori</i> Antibody	WB/S/P	Cassette	25 T	CE
IHP-601	<i>H. pylori</i> Antigen	Feces	Dipstick	25 T	CE
IHP-602	<i>H. pylori</i> Antigen	Feces	Cassette	25 T	CE
IVCA-602	<i>Vibrio cholerae</i> O1 (VC O1)	Feces	Cassette	10 T	CE
IVCB-602	<i>Vibrio cholerae</i> O139 (VC O139)	Feces	Cassette	10 T	CE
IVCAB-625	<i>Vibrio cholerae</i> O1/O139 Combo	Feces	Cassette	10 T	CE



Respiratory Tract Diseases

Cat No.	Product	Specimen	Format	Pack	CE Status
IST-501	Strep A (Control Line in Red)	Throat Swab	Dipstick	25 T	CE
IST-502	Strep A (Control Line in Red)	Throat Swab	Cassette	20 T	CE
ISTB-501	Strep A (Control Line in Blue)	Throat Swab	Dipstick	25 T	CE
ISTB-502	Strep A (Control Line in Blue)	Throat Swab	Cassette	20 T	CE
IAD-502	<i>Adenovirus pneumoniae</i>	Swab	Cassette	20 T	CE
IIN-501	Influenza A+B	Swab/Nasal Aspirate	Dipstick	25 T	CE
IIN-502	Influenza A+B	Swab/Nasal Aspirate	Cassette	20 T	CE
IFLA-501	Influenza A	Swab/Nasal Aspirate	Dipstick	25 T	CE
IFLA-502	Influenza A	Swab/Nasal Aspirate	Cassette	20 T	CE
IFA-502	H1N1 Antigen	Swab	Cassette	20 T	CE
IINA-525	Influenza A/B +H1N1 Combo	Swab	Cassette	10 T	CE
IRS-502	Respiratory Syncytial Virus Antigen (RSV)	Nasopharyngeal Swab /Nasal Aspirate	Cassette	20 T	CE
IRIC-525	RSV&Influenza A+B Combo	Swab/Nasal Aspirate	Cassette	10 T	CE
IARC-525	Adenovirus&RSV Combo	Nasopharyngeal Swab	Cassette	10 T	CE
IARI-535	Adenovirus, RSV and Influenza A+B Combo	Nasopharyngeal Swab	Cassette	10 T	CE
ILEG-102	<i>Legionella pneumophila</i>	Urine	Cassette	25 T	CE
ISP-102	<i>Streptococcus pneumoniae</i> Antigen	Urine	Cassette	10 T	CE
ISLC-125	<i>Streptococcus pneumoniae</i> and <i>Legionella pneumophila</i> Combo	Urine	Cassette	10 T	CE
IMP-502	<i>Mycoplasma pneumoniae</i> Antigen	Swab	Cassette	20 T	CE
IMPM-402	<i>Mycoplasma pneumoniae</i> IgM	WB/S/P	Cassette	25 T	CE
IMP-425	<i>Mycoplasma pneumoniae</i> IgG/IgM Combo	WB/S/P	Cassette	25 T	CE
IECG-402	EBV VCA IgG	WB/S/P	Cassette	10 T	CE
IENG-402	EBNA IgG	WB/S/P	Cassette	10 T	CE
IECNG-425	EBV VCA and EBNA IgG Combo	WB/S/P	Cassette	10 T	CE
IMO-402	MONO	WB/S/P	Cassette	25 T	CE
ICRG-402	<i>Cryptococcus</i> Antigen	WB/S/P/CSF	Cassette	10 T	CE



Rapid Tests for Infectious Diseases

Respiratory Tract Diseases - COVID-19

COVID-19 Antigen (For Professional)

Cat No.	Product	Specimen	Format	Pack	CE Status
INCP-ACO502	COVID-19 Antigen	Nasopharyngeal Swab	Cassette	20 T	CE
ISCO-ACO502	SARS-CoV-2 Antigen	Nasal Swab			
INCP-ACO802	COVID-19 Antigen	Oral fluid			

White List

Germany BFARM: AT930/21	Slovenia	France	Chile	Italy
Switzerland	Portugal	Austria	Turkey	Brazil


In addition, we have registered in more than 10 other countries, including Georgia, Ukraine, Indonesia, Argentina, Ecuador, Peru, etc.

COVID-19 Antigen (For Self-testing)

Cat No.	Product	Specimen	Format	Pack	CE Status
ACO-1001	SARS-CoV-2 Antigen	Nasal Swab	Cassette	1,2,5,7,10, 15,20,25 T	CE1434
ACO-1002	COVID-19 Antigen	Oral fluid			

Self-Testing List

CE1434	Germany BFARM: AT021/22	Switzerland	Austria
Brazil	Belgium	Lithuania	Portugal



• COVID-19 Antigen Rapid Test

- Easy specimen collection
- Fast results in 15 minutes assay time
- Easy visual interpretation
- Simple operation, no equipment required
- High accuracy

• COVID-19 Antigen Combo

Cat No.	Product	Specimen	Format	Pack	CE Status
ICIC-ACO525	COVID-19 and Influenza A+B Antigen Combo	Nasopharyngeal Swab	Cassette	20 T	CE
ISIR-ACO535	SARS-CoV-2/Influenza A+B/RSV Antigen Combo				
IRT-ACO545	SARS-CoV-2/Influenza A+B/RSV/Adenovirus Antigen Combo				
IRT-ACO555	SARS-CoV-2/Influenza A+B/RSV/Adenovirus/M.pneumoniae Antigen Combo				

• COVID-19 Antibody

Cat No.	Product	Specimen	Format	Pack	CE Status
INCP-402	2019-nCoV IgG/IgM	WB/S/P	Cassette	25 T	CE
INCP-402B	2019-nCoV IgG/IgM	WB/S/P	Cassette	25 T	CE
INCP-411B	COVID-19 IgG/IgM	WB/S/P	Dipstick	200 T	CE
INCPG-402	COVID-19 IgG	WB/S/P	Cassette	25 T	CE
INCPG-411	COVID-19 IgG	WB/S/P	Dipstick	200 T	CE
ICN-ACO402	COVID-19 Nabs	WB/S/P	Cassette	25 T	CE

• COVID-19 Antibody Combo

Cat No.	Product	Specimen	Format	Pack	CE Status
INMC-425	COVID-19 and M.pneumoniae IgG/IgM Combo	WB/S/P	Cassette	25 T	CE
INMCG-425	COVID-19 and M.pneumoniae IgG Combo				

Rapid Tests for Infectious Diseases


Hepatitis Markers

Cat No.	Product	Specimen	Format	Pack	CE Status
IHA-402	HAV IgM	WB/S/P	Cassette	25 T	CE
IHAGM-425	HAV IgG/IgM Combo	WB/S/P	Cassette	25 T	CE
IHBSG-301	HBsAg	S/P	Dipstick	50 T	Non-CE
IHBSG-302	HBsAg	S/P	Cassette	25 T	Non-CE
IHBSG-401	HBsAg	WB/S/P	Dipstick	50 T	Non-CE
IHBSG-402	HBsAg	WB/S/P	Cassette	25 T	Non-CE
IHBSB-401	HBsAb	WB/S/P	Dipstick	50 T	Non-CE
IHBSB-402	HBsAb	WB/S/P	Cassette	25 T	Non-CE
IHBSB-301	HBsAb	S/P	Dipstick	50 T	Non-CE
IHBSB-302	HBsAb	S/P	Cassette	25 T	Non-CE
IHBCB-302	HBcAb	S/P	Cassette	25 T	Non-CE
IHBEB-302	HBeAb	S/P	Cassette	25 T	Non-CE
IHBEG-302	HBeAg	S/P	Cassette	25 T	Non-CE
IHB-355	HBV Combo	S/P	Cassette	25 T	Non-CE
IHC-301	HCV	S/P	Dipstick	50 T	Non-CE
IHC-302	HCV	S/P	Cassette	25 T	Non-CE
IHC-401	HCV	WB/S/P	Dipstick	50 T	Non-CE
IHC-402	HCV	WB/S/P	Cassette	25 T	Non-CE
IHE-302	HEV IgG/IgM	S/P	Cassette	25 T	CE
IHBC-325	HBsAg and HCV Combo	S/P	Cassette	25 T	Non-CE
IHBC-425	HBsAg and HCV Combo	WB/S/P	Cassette	25 T	Non-CE



Blood Banking & Sexually Transmitted Infections (STIs)


Cat No.	Product	Specimen	Format	Pack	CE Status
IHI-302	HIV 1.2	S/P	Cassette	25 T	Non-CE
IHI-401	HIV 1.2	WB/S/P	Dipstick	50 T	Non-CE
IHI-402	HIV 1.2	WB/S/P	Cassette	25 T	Non-CE
IHI-402S	HIV 1.2 (Single Use Kit)	WB	Cassette	1 T	Non-CE
IHI-T401	HIV 1.2.O	WB/S/P	Dipstick	50 T	Non-CE
IHI-T402	HIV 1.2.O	WB/S/P	Cassette	25 T	Non-CE
IHI-T402S	HIV 1.2.O (Single Use Kit)	WB	Cassette	1 T	Non-CE
IHIG-402	HIV p24 Antigen	WB/S/P	Cassette	25 T	Non-CE
IHIBG-425	HIV 1.2 and HIV p24 Combo	WB/S/P	Cassette	25 T	Non-CE
IHIBG-402	HIV 1.2 and HIV p24 Combo (1 Window)	WB/S/P	Cassette	25 T	Non-CE
ISY-301	Syphilis	S/P	Dipstick	50 T	CE
ISY-302	Syphilis	S/P	Cassette	25 T	CE
ISY-401	Syphilis	WB/S/P	Dipstick	50 T	CE
ISY-402	Syphilis	WB/S/P	Cassette	25 T	CE
IMID-325	Syphilis/HIV1.2 Combo	S/P	Cassette	25 T	Non-CE
IISC-325	Syphilis/HIV1.2.O Combo	S/P	Cassette	25 T	Non-CE
IISC-425	Syphilis/HIV1.2.O Combo	WB/S/P	Cassette	25 T	Non-CE
IISC-402	HIV1.2 and Syphilis Combo (1 Window)	WB/S/P	Cassette	25 T	Non-CE
IBCH-335	HBsAg/HCV/HIV Combo	S/P	Cassette	25 T	Non-CE
IMID-345	HBsAg/HCV/HIV/Syphilis Combo	S/P	Cassette	25 T	Non-CE
IMID-445	HBsAg/HCV/HIV/Syphilis Combo	WB/S/P	Cassette	25 T	Non-CE
ITV-502	<i>Trichomonas Vaginalis</i>	Vaginal Swab	Cassette	10 T	CE
IGO-502	Gonorrhea	Swab	Cassette	20 T	CE
ICH-502	Chlamydia	Swab / Urine	Cassette	20 T	Non-CE
IGC-525	Gonorrhea and Chlamydia Combo	Swab	Cassette	10 T	Non-CE
ICA-502	<i>Candida albicans</i>	Swab	Cassette	10 T	CE
ISB-501	Strep B	Swab	Dipstick	25 T	CE
ISB-502	Strep B	Swab	Cassette	20 T	CE



Rapid Tests for Infectious Diseases


Tropical Diseases

Cat No.	Product	Specimen	Format	Pack	CE Status
ICHA-302	Chagas	S/P	Cassette	25 T	CE
ICHA-402	Chagas	WB/S/P	Cassette	25 T	CE
ICHI-302	Chikungunya IgG/IgM	S/P	Cassette	25 T	CE
ICHI-402	Chikungunya IgG/IgM	WB/S/P	Cassette	25 T	CE
IDEN-402	Dengue IgG/IgM	WB/S/P	Cassette	10 T	CE
IDES-402	Dengue NS1	WB/S/P	Cassette	10 T	CE
IDEC-425	Dengue IgG/IgM and NS1 Combo	WB/S/P	Cassette	10 T	CE
IZIG-402	Zika NS1	WB/S/P	Cassette	10 T	CE
IZIB-402	Zika IgG/IgM	WB/S/P	Cassette	10 T	CE
IZIC-425	Zika IgG/IgM and NS1 Combo	WB/S/P	Cassette	10 T	CE
IFIL-302	Filariasis IgG/IgM	S/P	Cassette	25 T	CE
IFIL-402	Filariasis IgG/IgM	WB/S/P	Cassette	25 T	CE
ITY-302	Typhoid IgG/IgM	S/P	Cassette	25 T	CE
ITY-402	Typhoid IgG/IgM	WB/S/P	Cassette	25 T	CE
ITY-401	Typhoid IgG/IgM	WB/S/P	Dipstick	50 T	CE
ISTY-602	<i>Salmonella Typhi</i> Antigen	Feces	Cassette	25 T	CE
ISPT-602	<i>Salmonella paratyphi</i> Antigen	Feces	Cassette	25 T	CE
ISTP-625	<i>Salmonella typhi</i> and <i>paratyphi</i> Antigen Combo	Feces	Cassette	10 T	CE
ILEI-402	Leishmania IgG/IgM	WB/S/P	Cassette	25 T	CE
ILEP-402	Leptospira IgG/IgM	WB/S/P	Cassette	25 T	CE
IMA-401	Malaria P.f.	WB	Dipstick	50 T	CE
IMA-402	Malaria P.f.	WB	Cassette	25 T	CE
IMPV-402	Malaria P.f./P.v.	WB	Cassette	25 T	CE
IMPV-402	Malaria P.f./P.v.	WB	Cassette	25 T	CE
IMPVF-402	Malaria P.f./P.v./Pan	WB	Cassette	25 T	CE




Others

Cat No.	Product	Specimen	Format	Pack	CE Status
IBA-402	Brucella Abortus Antigen	WB/S/P	Cassette	25 T	CE
ILY-402	Lyme IgG/IgM	WB/S/P	Cassette	10 T	CE
ITSU-402	Scrub Typhus IgG/IgM	WB/S/P	Cassette	10 T	CE
ITE-402	Tetanus	WB/S/P	Cassette	25 T	CE
ITB-401	Tuberculosis (TB)	WB/S/P	Dipstick	50 T	CE
ITB-402	Tuberculosis (TB)	WB/S/P	Cassette	25 T	CE
IMXB-ACO402	Monkeypox IgG/IgM	WB/S/P	Cassette	25 T	Non-CE



Rapid Tests for Drugs of Abuse

Test with Urine Specimen




Analyte	Cut-off	CE Status
AB-PINACA (ABP)	10 ng/mL	CE
Acetaminophen (ACE)	5000 ng/mL	CE
Alcohol (ALC)	0.02%	CE
7-Aminoclonazepam (7-ACL)	300 ng/mL	CE
Alprazolam (ALP)	100 ng/mL	CE
Amphetamine (AMP)	300 ng/mL, 500 ng/mL, 1000 ng/mL*	CE/US FDA
Δ-Pyrrolidinovalerophenone (Δ-PVP)	1000 ng/mL*, 500 ng/mL, 300 ng/mL	CE
Barbiturate (BAR)	200 ng/mL, 300 ng/mL*	CE/US FDA
Buprenorphine (BUP)	5 ng/mL, 10 ng/mL*	CE/US FDA
Benzodiazepines (BZO)	100 ng/mL, 200 ng/mL, 300 ng/mL*	CE/US FDA
Cathine (CAT)	150 ng/mL	CE
Caffeine (CAF)	1000 ng/mL	CE
Carisoprodol (CAR)	2000 ng/mL	CE
Cannabinol (CNB)	500 ng/mL	CE
Carfentanyl (CFYL)	500 ng/mL*, 250 ng/mL	CE
Citalopram (CIT)	500 ng/mL	CE
Clonazepam (CLO)	150 ng/mL, 400 ng/mL*	CE
Cocaine (COC)	100 ng/mL, 150 ng/mL, 300 ng/mL*	CE/US FDA
Cotinine (COT)	10 ng/mL, 50 ng/mL, 100 ng/mL, 200 ng/mL*, 500 ng/mL	CE
Diazepam (DIA)	300 ng/mL	CE
Ethylendiamine-dimethylphosphinic acid (EDDP)	300 ng/mL*, 100 ng/mL	CE
Ethyl Glucuronide (ETG)	500 ng/mL*, 1000 ng/mL	CE
F-Ketamine (FKET)	1000 ng/mL*, 500 ng/mL	CE
Fluoxetine (FLX)	500 ng/mL	CE
Fentanyl (FYL)	10 ng/mL, 20 ng/mL*, 100 ng/mL, 200 ng/mL	CE
Gabapentin (GAB)	2000 ng/mL	CE
γ-Hydroxybutyric acid (GHB)	10 µg/mL	CE
Ketamine (KET)	100 ng/mL, 300 ng/mL, 500 ng/mL, 1000 ng/mL*	CE
Kratom (KRA)	300 ng/mL	CE
Lysergic Acid Diethylamide (LSD)	10 ng/mL, 20 ng/mL*, 50 ng/mL	CE
6-Monoacetylmorphine (6-MAM)	10 ng/mL	CE
Methcathinone (MCAT)	500 ng/mL	CE
Ecstasy (MDMA)	500 ng/mL*, 1000 ng/mL	CE/US FDA
Tenamfetamine (MDA)	500 ng/mL	CE
Methylenedioxypyrovalerone (MDPV)	1000 ng/mL*, 500ng/mL, 300ng/mL	CE

Analyte	Cut-off	CE Status
Methamphetamine (MET)	300 ng/mL, 500 ng/mL, 1000 ng/mL*	CE/US FDA
Mephedrone HCl (MEP)	100 ng/mL	CE
Mescaline (MES)	100 ng/mL*, 300ng/mL	CE
Morphine (MOP)	100 ng/mL, 300 ng/mL*	CE/US FDA
Methylphenidate (MPD)	1000 ng/mL*, 300 ng/mL	CE
Meperidine (MPRD)	100 ng/mL	CE
Methaqualone (MQL)	300 ng/mL	CE
Methadone (MTD)	100 ng/mL, 200 ng/mL, 300 ng/mL*	CE/US FDA
Mirtazapine (MZT)	500 ng/mL	CE
N, N-Dimethyltryptamine (NND)	1000 ng/mL	CE
Opiates (OPI)	2000 ng/mL*, 300 ng/mL	CE/US FDA
Oxycodone (OXY)	100 ng/mL	CE/US FDA
Olanzapine (OZP)	1000 ng/mL	CE
Papaverine (PAP)	500 ng/mL	CE
Phencyclidine (PCP)	25 ng/mL	CE/US FDA
Pregabalin (PGB)	50 µg/mL*, 500 ng/mL	CE
Propoxyphene (PPX)	300 ng/mL	CE
Quetiapine (QTP)	1000 ng/mL	CE
Risperidone (RPD)	150 ng/mL	CE
Scopolamine (SCOP)	500 ng/mL	CE
Synthetic Marijuana (K2)	50 ng/mL*, 30 ng/mL	CE
Tapentadol (TAP)	1000 ng/mL	CE
Tricyclic Antidepressants (TCA)	500 ng/mL, 1000 ng/mL*	CE/US FDA
Marijuana (THC)	20 ng/mL, 25 ng/mL, 50 ng/mL*, 150 ng/mL, 300 ng/mL	CE/US FDA
Tilidine (TLD)	50 ng/mL	CE
Tramadol (TML)	100 ng/mL*, 200 ng/mL, 300 ng/mL	CE
Tropicamide (TRO)	350 ng/mL	CE
Trazodone (TZD)	200 ng/mL	CE
UR-144	25 ng/mL	CE
Zaleplon (ZAL)	100 ng/mL	CE
Zolpidem (ZOL)	50 ng/mL	CE
Zopiclone (ZOP)	50 ng/mL	CE


* Standard Cut-off

Formats & Pack : Dipstick (50 T), Cassette (25 T), Panel (25 T), Cup (25 T)

Multi-Drug 2~17 Drugs Rapid Test Cassette



Multi-Drug 2~20 Drugs Rapid Test Panel



Multi-Drug 2~22 Drugs Rapid Test 1-Step Cup



Multi-Drug 2~16 Drugs Rapid Test Key Cup

Rapid Tests for Drugs of Abuse

Test with Oral Fluid Specimen

Analyte	Cut-off	CE Status
Amphetamine (AMP)	50 ng/mL	CE
Δ-Pyrrolidinovalerophenone (Δ-PVP)	300 ng/mL	CE
Barbiturates (BAR)	50 ng/mL	CE
Buprenorphine (BUP)	10 ng/mL	CE
Benzodiazepines (BZO)	20 ng/mL, 30 ng/mL*, 50 ng/mL	CE
Cocaine (COC)	20 ng/mL	CE
Carfentanyl (CFYL)	50 ng/mL	CE
Cotinine (COT)	20 ng/mL	CE
Fentanyl (FYL)	50 ng/mL*, 20 ng/mL	CE
Ketamine (KET)	50 ng/mL	CE
6-Monoacetylmorphine (6-MAM)	10 ng/mL	CE
Ecstasy (MDMA)	50 ng/mL	CE
Methylenedioxypyrovalerone (MDPV)	300 ng/mL	CE
Methamphetamine (MET)	50 ng/mL	CE
Methadone (MTD)	30 ng/mL	CE
Opiates (OPI)	40 ng/mL	CE
Oxycodone (OXY)	20 ng/mL	CE
Phencyclidine (PCP)	10 ng/mL	CE
Synthetic Marijuana (K2)	25 ng/mL	CE
Marijuana (THC) (Parent)	20 ng/mL, 50 ng/mL*	CE
Marijuana (THC)	50 ng/mL	CE
Tramadol (TML)	30 ng/mL	CE

* Standard Cut-off

Formats & Pack : Cassette (25 T), Device (25 T), Cup (20 T)



Multi-Drug 2-12 Drugs Rapid Test Device



Multi-Drug 2-9 Drugs Rapid Test Cassette




Multi-Drug 2-16 Drugs
Rapid Test C Cup

Test with Hair Specimen

Analyte	Cut-off	CE Status
AB-PINACA (ABP)	0.2 ng/mg	CE
Amphetamine (AMP)	0.2 ng/mg	CE
Barbiturate (BAR)	0.2 ng/mg	CE
Buprenorphine (BUP)	0.2 ng/mg	CE
Benzodiazepines (BZO)	0.2 ng/mg	CE
Cocaine (COC)	0.5 ng/mg	CE
Cotinine (COT)	0.2 ng/mg	CE
Ketamine (KET)	0.2 ng/mg	CE
6-Monoacetylmorphine (6-MAM)	0.2 ng/mg	CE
Ecstasy (MDMA)	0.2 ng/mg	CE
Methamphetamine (MET)	0.2 ng/mg	CE
Morphine (MOP)	0.2 ng/mg	CE
Oxycodone (OXY)	0.2 ng/mg	CE
Phencyclidine (PCP)	0.2 ng/mg	CE
Synthetic Marijuana (K2)	0.2 ng/mg	CE

Formats & Pack : Cassette (10 T)



Rapid Tests for Drugs of Abuse

CE
Multi-Drug 2~16
Drugs Rapid Test
Cassette



Test with Powder Specimen (Forensic Use)

Analyte	Cut-off	CE Status
Acetaminophen (ACE)	5000 ng/mL	Non-CE
Amphetamine (AMP)	1000 ng/mL	Non-CE
Δ-Pyrrolidinovalerophenone (Δ-PVP)	1000 ng/mL	Non-CE
Barbiturate (BAR)	300 ng/mL	Non-CE
Buprenorphine (BUP)	10 ng/mL	Non-CE
Benzodiazepines (BZO)	300 ng/mL	Non-CE
Cocaine (COC)	300 ng/mL	Non-CE
Carfentanyl (CFYL)	500 ng/mL	Non-CE
Cotinine (COT)	200 ng/mL	Non-CE
Diazepam (DIA)	300 ng/mL	Non-CE
Fentanyl (FYL)	200 ng/mL	Non-CE
Ketamine (KET)	1000 ng/mL	Non-CE
Lysergic Acid Diethylamide (LSD)	20 ng/mL	Non-CE
Ecstasy (MDMA)	500 ng/mL	Non-CE
Tenamfetamine (MDA)	500 ng/mL	Non-CE
Methylenedioxypyrovalerone (MDPV)	1000 ng/mL	Non-CE
Methylphenidate (MPD)	300 ng/mL	Non-CE
Methamphetamine (MET)	1000 ng/mL	Non-CE
Morphine (MOP)	300 ng/mL	Non-CE
Methaqualone (MQL)	300 ng/mL	Non-CE
Methadone (MTD)	300 ng/mL	Non-CE
Oxycodone (OXY)	100 ng/mL	Non-CE
Opiate (OPI)	2000 ng/mL	Non-CE
Phencyclidine (PCP)	25 ng/mL	Non-CE
Propoxyphene (PPX)	300 ng/mL	Non-CE
Synthetic Marijuana (K2)	50 ng/mL	Non-CE
Tricyclic Antidepressants (TCA)	1000 ng/mL	Non-CE
Marijuana (THC)	50 ng/mL	Non-CE
Tramadol (TML)	100 ng/mL	Non-CE
Zolpidem(ZOL)	50 ng/mL	Non-CE


Formats & Pack : Panel (25 T)

Test with Whole Blood / Serum / Plasma Specimen

Analyte	Cut-off	CE Status
Amphetamine (AMP)	80 ng/mL	CE
Barbiturate (BAR)	100 ng/mL	CE
Buprenorphine (BUP)	5 ng/mL	CE
Benzodiazepines (BZO)	100 ng/mL	CE
Cocaine (COC)	50 ng/mL	CE
Cotinine (COT)	10 ng/mL, 100 ng/mL*	CE
Ethylendiamine-dimethylphosphinic acid (EDDP)	50 ng/mL	CE
Fentanyl (FYL)	15 ng/mL	CE
Ketamine (KET)	200 ng/mL	CE
Lysergic Acid Diethylamide (LSD)	20 ng/mL	CE
Ecstasy (MDMA)	50 ng/mL	CE
Tenamfetamine (MDA)	80 ng/mL	CE
Methylenedioxypyrovalerone (MDPV)	300 ng/mL	CE
Methamphetamine (MET)	70 ng/mL	CE
Morphine (MOP)	40 ng/mL	CE
Methadone (MTD)	40 ng/mL	CE
Oxycodone (OXY)	20 ng/mL	CE
Phencyclidine (PCP)	20 ng/mL	CE
Propoxyphene (PPX)	100 ng/mL	CE
Synthetic Marijuana (K2)	100 ng/mL	CE
Tricyclic Antidepressants (TCA)	300 ng/mL	CE
Marijuana (THC)	35 ng/mL*, 50 ng/mL	CE
Tramadol (TML)	50 ng/mL	CE


* Standard Cut-off

Formats & Pack : Panel (25 T)



Multi-Drug 2~17 Drugs Rapid Test Cassette

Rapid Tests for Tumor Markers



Cat No.	Product	Specimen	Format	Pack	Cut-off	CE Status
TAF-402	Alpha-Fetoprotein (AFP)	WB/S/P	Cassette	25 T	10 ng/mL	CE
TCE-402	Carcinoembryonic Antigen (CEA)	WB/S/P	Cassette	25 T	5 ng/mL	CE
TPS-Q401	Prostate Specific Antigen (PSA) (Qualitative)	WB/S/P	Dipstick	25 T	4 ng/mL	Non-CE
TPS-Q402	Prostate Specific Antigen (PSA) (Qualitative)	WB/S/P	Cassette	25 T	4 ng/mL	Non-CE
TPS-402	Prostate Specific Antigen (PSA)	WB/S/P	Cassette	25 T	3-10 ng/mL	Non-CE
T125-402	CA125	WB/S/P	Cassette	10 T	40 U/mL	CE
T153-402	CA15-3	WB/S/P	Cassette	10 T	30 U/mL	CE
T199-402	CA19-9	WB/S/P	Cassette	10 T	40 U/mL	CE
TFO-602	Fecal Occult Blood (FOB)	Feces	Cassette	25 T	25 ng/mL, 40 ng/mL, 50 ng/mL*, 100 ng/mL	CE
TFO-601	Fecal Occult Blood (FOB)	Feces	Dipstick	25 T		CE
TCFC-625	Calprotectin and FOB Combo	Feces	Cassette	10 T	See Insert	CE
TTFC-625	Transferrin and FOB Combo	Feces	Cassette	25 T	See Insert	CE
THPB-625	Hb+Hb-Hp Combo	Feces	Cassette	25 T	50 ng/mL	CE
TTFH-635	Transferrin/FOB and Hb-Hp Combo	Feces	Cassette	25 T	See Insert	CE

Cardiac Markers

Cat No.	Product	Specimen	Format	Pack	Cut-off	CE Status
CTI-402	Cardiac Troponin I (cTnI)	WB/S/P	Cassette	25 T	0.5 ng/mL	CE
CTNT-402	Cardiac Troponin T (cTnT)	WB/S/P	Cassette	10 T	0.5 ng/mL	CE
CCK-402	CK-MB	WB/S/P	Cassette	25 T	5 ng/mL	CE
CMY-402	Myoglobin	WB/S/P	Cassette	25 T	50 ng/mL	CE
CFA-402	H-FABP	WB/S/P	Cassette	10 T	8 ng/mL	CE
CFC-425	H-FABP and cTnI Combo	WB/S/P	Cassette	10 T	See Insert	CE
CMA-435	Myoglobin/CK-MB/Troponin I Combo	WB/S/P	Cassette	25 T	See Insert	CE
CMA-445	H-FABP and Myo/CK-MB/cTnI Combo	WB/S/P	Cassette	10 T	See Insert	CE
CBN-402	NT-proBNP	WB/S/P	Cassette	10 T	0.45 U/mL	CE
CDM-402	D-dimer	WB/S/P	Cassette	10 T	500 ng/mL	CE
CCR-401	C-reactive protein (CRP)	WB/S/P	Dipstick	20 T	10 µg/mL	CE
CCR-402	C-reactive protein (CRP)	WB/S/P	Cassette	20 T	10 µg/mL*, 0.5 µg/mL	CE
CCR-U401	C-reactive protein (CRP) (Semi-Quantitative)	WB/S/P	Dipstick	20 T	10-40-80 µg/mL	CE
CCR-T402	C-reactive protein (CRP) (Semi-Quantitative)	WB/S/P	Cassette	20 T	10-30 µg/mL	CE
CCR-U402	C-reactive protein (CRP) (Semi-Quantitative)	WB/S/P	Cassette	20 T	10-40-80 µg/mL	CE
CPC-302	Procalcitonin (PCT)	S/P	Cassette	10 T	1 ng/mL	CE
CPC-402	Procalcitonin (PCT)	WB/S/P	Cassette	10 T	1 ng/mL	CE

Others

Cat. No.	Product	Specimen	Format	Pack	Cut-off	CE Status
OFE-402	Ferritin	WB/S/P	Cassette	10 T	30 ng/mL	CE
OFE-T402	Ferritin (Semi-Quantitative)	WB/S/P	Cassette	10 T	13-30 ng/mL	CE
OSP-902	SP-10 Male Fertility	Sperm	Cassette	1 T/2 T	15 million/mL	CE
OTS-402	TSH	WB/S/P	Cassette	25 T	5 μ IU/mL	CE
OVD-402	Vitamin D	WB	Cassette	10 T	10-30-100 ng/mL	CE
OHA-402	HbA1c	WB	Cassette	10 T	6%~9%	CE
OBS-901	Blood Stain	Stain	Dipstick	1 T	50 ng/mL	Non-CE
OBS-902	Blood Stain	Stain	Cassette	10 T	50 ng/mL	Non-CE
OHS-901	Human Semen	Stain	Dipstick	50 T	4 ng/mL	Non-CE
OHS-902	Human Semen	Stain	Cassette	25 T	4 ng/mL	Non-CE
OCAL-602	Calprotectin	Feces	Cassette	10 T	50 μ g/g	CE
OLF-602	Lactoferrin	Feces	Cassette	10 T	100 ng/mL	CE
OIBD-625	Calprotectin and Lactoferrin Combo	Feces	Cassette	10 T	/	CE
OIGE-402	IgE	WB/S/P	Cassette	10 T	200 IU/mL	CE
ODME-ACO402	Dust Mite IgE	WB/S/P	Cassette	10 T	0.35 IU/mL	CE
ORF-302	Rheumatoid Factor	S/P	Cassette	20 T	25 IU/mL	CE
OMIA-101	Micro-Albumin (Semi-Quantitative)	Urine	Dipstick	50 T	20-40-60-80 μ g/mL	CE
OMIA-102	Micro-Albumin (Semi-Quantitative)	Urine	Cassette	25 T	20-40-60-80 μ g/mL	CE
OMAL-101	Micro-Albumin (Qualitative)	Urine	Dipstick	50 T	20 μ g/mL	CE
OMAL-102	Micro-Albumin (Qualitative)	Urine	Cassette	25 T	20 μ g/mL	CE
OSA-402	SAA	WB/S/P	Cassette	10 T	10 μ g/mL	CE
OCSC-425	SAA and CRP Combo	WB/S/P	Cassette	10 T	10 μ g/mL	CE
OAB-402	ABO Blood Grouping	WB	Cassette	25 T	/	Non-CE
OABD-402	ABO and RhD Blood Grouping	WB	Cassette	25 T	/	Non-CE



NOTE:

NOTE:

The logo for JusChek+ features the word "JusChek" in a bold, black, sans-serif font, with a small orange plus sign superscripted to the right of the "k".

JusChek⁺

Version:11



Via XXI Ottobre 1944, 11/2 - 40055 Castenaso (BO)
Tel.: 051.6240700 - Fax: 051.6240706
e-mail: diagnostica@temaricerca.com
web: www.temaricerca.com