

# Fluorescent Immunoanalyzer

## AFS-1000(Single Channel)

One-step operation to complete the testing. Automatically identify test items, lost card, print result, and upload to LIS.

### Application:

Various medical institutions(clinical departments, emergency departments , community medical centers , primary health centers). Vet hospital, food safety inspection, inspection and quarantine, environmental monitoring , etc.

### Certification:

CE &FDA



### //Technical performance//

• Principle	Dry fluorescence immunoassay
• Channel	Single
• Sample Type	Whole blood/Serum/Plasma/Urine/Others
• Reading Time	<10s
• Testing Method	Support one reagent card with multiple test items
• Reagent Card Code	Bar code& QR code
• Printer	Embedded thermal printer
• Reagent	Customized: Conventional card, Microfluidic chip card
• Transmission	4G, USB2.0, RS232C, WIFI, Ethernet, LIS, HIS
• Screen	7 inches Touch screen
• Language	English and others
• Size	215mm*302mm*155mm
• Weight	3Kg
• Data Storage	Unlimited